

# RABTA

Random Access Biochip Technology Analyser



# INNOVATION MEETS PRECISION



# CONTENTS

Biochip Array Technology

Why RABTA?

RABTA Modules

Analyser Workflow

Test Menu

Software Overview

Technical Snapshot

# **RABTA**

The new RABTA (Random Access Biochip Technology Analyser) sets a high standard in laboratory diagnostics by combining advanced biochip technology with exceptional efficiency. Designed for seamless operation, the RABTA platform allows uninterrupted workflows and increased productivity. With rapid random-access capabilities, the analyser provides precise and reliable results across a wide range of assays.




Featuring five intuitive modules, this innovative analyser simplifies high-throughput biological analysis, enabling clinicians to conduct complex assays quickly and accurately. Each module is tailored for specific functions, ensuring easy integration and adaptability for various applications, from diagnostics to advanced research. Its random-access capability allows simultaneous processing of multiple samples, making it an essential tool for modern laboratories.

# The Power of Biochip Technology

Randox Biochip Array Technology is an innovative multiplex testing platform that allows for the simultaneous detection of multiple analytes from a single sample. This technology is based on a remarkable 9x9mm biochip, which functions as a solid-phase reaction vessel. Each biochip is precoated with an array of discrete test regions (DTRs), enabling the detection of up to 48 different tests from one sample.

Using competitive chemiluminescent immunoassays when a sample is applied to the biochip, specific reactions occur at each DTR, producing light signals. These signals are detected and analysed using sophisticated digital imaging technology revealing the concentration of each analyte with pinpoint accuracy.



By consolidating multiple assays onto one biochip, the Biochip Array reduces the time and resources needed for diagnostics, making it an indispensable tool for modern diagnostics. This efficiency has transformed how laboratories are analysing samples, processing a higher volume without compromising speed or accuracy.

# **Applications**



Research Laboratories



Diagnostics Laboratories



Public & Private Hospitals



University Laboratories



Private Laboratories



Commercial Laboratories



Forensic/Criminal Laboratories



Clinical Laboratories



# WHY RABTA?



#### High Throughput

A maximum throughput of up to 60 samples can be loaded per hour, providing up to 2,640 test results per hour / 19,580 test results per day.



#### **Priority Sampling**

Software allows users to assign samples as a priority over others so that **priority samples** can be run and results reported rapidly, **eliminating the constraints of batch testing.** 



#### Workflow Consolidation

The RABTA consolidates the workflow of multiple laboratories into a single analyser. The workload of laboratory staff is minimised, as only consumables and samples need to be loaded, and tests are selected via onscreen software.



#### 2.5 Hour Walk Away

A fully automated platform with an operator walkaway time of up to 2.5 hours. A single operator is all that's required to run the RABTA.



#### Zero Carry-Over Risk

With aspiration abilities from 1.5µl upward, the RABTA does not need pre-dilution. The system utilises single-use tips for aspiration and dispense, meaning there is no carry-over risk at all, ensuring high accuracy with every result.



#### **Module Access**

Continuous access to samples.
Software-controlled access to reagents, biochips, signal and waste.
The system can hold up to 60 samples with continuous access, and 120 biochips with software-controlled access.



#### **Multiplex Testing**

Simultaneous detection of up to 44 targets from a single patient sample through our **patented Biochip Array Technology.** 



#### **Rapid Results**

Software intuitively analyses what tests are being run and assigns the order of samples in the most time-efficient way for analysing and reporting results. The time to first sample results is just 36 minutes, with up to 44 results per sample. Results for remaining samples are released 1 sample per minute onward.



#### Rapid Start Up

The system will be ready for use in under 15 minutes when started from a standby state.



# **RABTA MODULES**



# Reagent Storage

#### **KEEPS REAGENTS FRESH AND READY**

Onboard storage chills reagents at a perfect 2-8°C, ensuring they stay stable and effective for up to 7 days. No more worries about reagent degradation.



# **Biochip Storage**

#### **MAXIMISE EFFICIENCY**

Holding up to 120 biochips through softwarecontrolled access, the RABTA ensures uninterrupted workflows and improves everyday efficiency.



# Signal Storage

#### **POWER UP ASSAYS**

Peroxide and luminol are securely stored in this module, ready to deliver the signals needed for clear, reliable results. Assays are analysed with seamless precision and consistency.



## **Tip Storage**

#### **NO CROSS-CONTAMINATION**

The RABTA platform can store up to 768 single-use tips, each one dedicated to precise aspiration and dispensing of assay diluent, samples and conjugates onto biochips. This means every result is accurate, every time.



### Sample Storage

#### **FLEXIBILITY AT YOUR FINGERTIPS**

With rapid loading and access, samples can be added for quick and easy analysis. Load up to 60 samples and let the platform handle the rest, ensuring seamless operation and timely reporting.

# **Optimised Waste Management**

#### **Custom Bench**

The RABTA platform includes a specially designed Custom Bench. This not only provides ample storage space for all your consumables but also ensures smooth and efficient disposal of both solid and liquid waste.



#### **Optimised Waste Management**

Featuring advanced sensors within the custom bench the waste bins alert the user when it's time to empty. With a generous waste storage capacity of up to 20 litres, there are longer intervals between emptying, significantly reducing hands-on time and streamlining the user's workflow.

# **ANALYSER WORKFLOW**

System auto primes in under 10 mins.



Consumables loaded, e.g. tips, biochips, signal and reagents



Biochip peeled and placed into dispense seat



Samples loaded



ADP aspirates and dispenses assay diluent and sample into a well



User loads tests in the software



Biochip is placed into the incubator



Press start



Conjugate dispensed into the well whilst inside the incubator





Incubated for set time



System takes and processes images



Moved to wash station and washed followed by spin dry



Biochip moved to waste



Moved to signal seat and 500µl signal dispensed



Results generated on screen and exported to LIMS



Moved to image station

# TEST MENU



Male & Female Hormones



Lifestyle & Nutrition



Toxicology



Oncology & Immunology



Haematology & Cardiology

RABTA 0

#### Male & Female Hormones

#### Fertility Biochip







Serum



Sample Volume:  $125 \mu l$ 



Time to Result:

36 minutes
(first sample result)



Result Type: Quantitative

# Biomarkers Prolactin

Prolactin Luteinizing Hormone (LH) Follicle Stimulating Hormone (FSH) Estradiol Progesterone

#### SHBG Biochip |



Sample Type: Serum



Sample Volume:  $10 \mu l$ 



Time to Result:





Result Type:

Quantitative



SHBG (Sex Hormone-Binding Globulin)

#### **Testosterone Biochip**



Sample Type: Serum



Sample Volume:  $125 \mu l$ 



Time to Result:

36 minutes (first sample result)



#### Result Type:

Quantitative

#### Biomarkers

Total Testosterone

#### Lifestyle & Nutrition

#### Vitamin D Biochip



Sample Type: Serum



Sample Volume:  $90 \mu l$ 



Time to Result:

36 minutes (first sample result)



Result Type:

Quantitative



25-hydroxy vitamin D (total)

#### **Toxicology**

#### ToxPlex Urine





Sample Type:

Urine



Time to Result:

36 minutes (first sample result)



Sample Volume:

25 μΙ



Result Type:

Semi-Quantitative

#### **Biomarkers**

Methamphetamine MDMA Amphetamine Tricyclic Antidepressants Oxycodone

Opiate 6-Monoacetylmorphine (6-MAM) Benzodiazepines 1 Benzodiazepines 2

Barbiturates

Cannabinoids (THC) Dextromethorphan Methadone Benzoylecgonine (BZG) / Cocaine Meprobamate Tramadol

Fentanvl Buprenorphine Propoxyphene Phencyclidine (PCP) Zolpidem

Ketamine Haloperidol Methaqualone Pregabalin Creatinine Ethyl Glucuronide (EtG) Acetaminophen Salicylates

#### Drugs of Abuse VII Urine





Sample Type:

Urine





Time to Result:

36 minutes (first sample result) Sample Volume:

25 μΙ



Result Type:

Semi-Quantitative

#### Biomarkers

Amphetamine Barbiturates Benzoylecgonine (BZG) / Benzodiazepines Opiate Cannabinoids (THC) Methamphetamine

Tramadol Creatinine

#### Drugs of Abuse VI Urine



17



Sample Type:

Urine





Time to Result:

36 minutes (first sample result)



Sample Volume:

25 μΙ



Result Type:

Semi-Quantitative

**Biomarkers** 

AB-PINACA Amphetamine Barbiturates Benzodiazepines Buprenorphine Creatinine

Benzoylecgonine (BZG) / Cocaine Cannabinoids (THC) Ethyl Glucuronide (EtG) JWH018

MDMA Methamphetamine Opiates UR144 6-Monoacetylmorphine (6-MAM)

Research Use Only

#### Oncology & Immunology

#### Carcinoembryonic Antigen (CEA) Biochip I





Sample Type:

Serum





Sample Volume: 100 μΙ



Time to Result: 36 minutes (first sample result)



Result Type: Quantitative

Biomarkers

Carcinoembryonic Antigen (CEA)

#### PD-1 Biochip



Sample Type: Serum



Sample Volume: 100 μΙ



Time to Result: 36 minutes (first sample result)



Result Type: Quantitative



Programmed Cell Death Protein 1

#### Mesothelin Biochip



Sample Type: Serum



Sample Volume: 100 μΙ



Time to Result: 36 minutes (first sample result)



Result Type: Quantitative

Biomarkers Mesothelin

Research Use Only 18

#### Alpha-fetoprotein (AFP) Biochip I











Biomarkers
Alpha-fetoprotein (AFP)

Sample Type: Serum Sample Volume:  $50 \mu l$ 

Time to Result:

36 minutes
(first sample result)

**Result Type:**Quantitative

## Haematology & Cardiology

#### Galectin-3 Biochip











Biomarkers Galectin-3

Sample Type: Sam
Serum


iample Volume: 20 μl Time to Result:

36 minutes
(first sample result)

**Result Type:**Quantitative


Research Use Only

# **SOFTWARE OVERVIEW**




#### Consumable Overview Data

Cutting-edge software offers a comprehensive view of all consumables onboard, ensuring users always know the sample status at a glance. Each section provides detailed insights for continuous monitoring that helps streamline workflows and boost laboratory efficiency. The system smartly prevents test processing if consumables are insufficient, prompting users to reload as needed. Plus, it halts tests if any consumables have expired, ensuring only the freshest materials are used.



#### **Analysis Overview**

The RABTA advanced software offers a comprehensive analysis overview, displaying all scheduled and running tests. Track progress effortlessly on the right-hand side, where users can see exactly how much time is left for each test. The software intelligently analyses ongoing tests and optimises the order of samples for maximum efficiency in analysis and reporting. The system allows you to assign priority to STAT samples, ensuring rapid processing and reporting of critical results.



#### Loaded Biochips Data

The dynamic onscreen cartridge overview displays all loaded biochips in real time, showing live biochip counts. Users can easily load, unload, and manually enter cartridge details with the intuitive interface. Cartridges are colour-coded and marked with their positions for effortless identification, making lab operations smoother and more efficient than ever before.

# Technical Snapshot

CAT Number EV4508

**Dimensions** 1300 x 800 x 1800mm (W x D x H)

Weight 450kg

Analyser Type Fully Automated Immunoassay platform

Biochips per Cartridge 10

Maximum Sample Throughput 60 samples per hour

Maximum Test Result Throughput 2,640 test results per hour

**Data Storage Capacity**A minimum of 1 year's worth of data can be stored

Walk Away Time 2.5 hours

Maximum Sample/Biochip Capacity

The system can hold up to 60 samples (6x10) with continuous access, and 100 biochips (12, 10) with a few years and to live the second samples (6x10) with a few years and to live the second samples (6x10) with a few years and to live the second samples (6x10) with a few years and to live the second samples (6x10) with a few years and to live the second samples (6x10) with continuous access, and the second samples (6x10) with a second sample (6x10) with a s

and 120 biochips (12x10) with software-controlled access.

Sample Types Serum, Urine

Time to First Result 36 minutes, 1 sample per minute onwards

Sample Clot Detection Yes

Priority Samples Yes, software controlled

Maximum Reagent/Assay Capacity on Board The system can hold up to 12 reagent types to match biochips

Continuous Access To Samples

Random Access Yes

Software-Controlled Access To Reagents, Biochips, Signal and Waste

QC Master Curve adjuster calibration along with control runs

**Calibration Time** 2 sample adjusters. 38 minutes plus processing time.

Calibration Frequency Required every 28 days

LIMS Compatible Uni-directional LIMs

Operational Modes Windows 10

Operating Temperature Room temp +18°C - +32°C

Daily – 5 mins

**Operator Maintenance** Weekly – 15mins

Monthly - 1 hour

Onboard Maintenance Records Yes

Remote Diagnostics Yes

**Emergency Stop Feature** Yes

Automatic Shutdown Feature Yes

Alerts Yes, alerts are shown onscreen





