



STI-*flex* qPCR Assays

qPCR assays available as four separate kits for increased flexibility as per laboratory requirements. Universal cycling conditions with PCR run-time <65 minutes to result. Internal Control for monitoring of sample and assay integrity throughout the process.

FEATURES



qPCR assays for the sensitive detection of 9 pathogens



Each kit has a 96-sample testing capacity



Rapid turnaround time with PCR run time <65 minutes



Validated for a range of sample matrices including urine and urogenital swabs

PATHOGENS DETECTED:

Chlamydia trachomatis (CT)

Neisseria gonorrhoeae (NG)

Herpes simplex virus 1 (HSV-1)

Herpes simplex virus 2 (HSV-2)

Treponema pallidum (*Syphilis*) (TP)

Mycoplasma genitalium (MG)

Trichomonas vaginalis (TV)

Mycoplasma hominis (MH)

Ureaplasma urealyticum (UU)

marketing@randox.com | www.randox.com



Copyright © 2025 Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151682708. Product availability may vary from country to country. Some products may be for Research Use Only. For more information on product application and availability, please contact your local Randox Representative

LT 1000 JUL25



BENEFITS

- » STI-*flex* qPCR Assays are offered in four separate kits to provide increased *flexibility*.
- » Design includes an Internal Control to prevent false negatives and ensure successful DNA extraction and PCR.
- » Efficient PCR set up with only 2 components required for mastermix.
- » Each STI-*flex* assay uses the same PCR cycling conditions; therefore, multiple assays can be run on the same instrument concurrently.

| Catalogue Code | Assay | Pack Size |
|----------------|-------------------------------------|-----------|
| QP10704 | CT/NG <i>flex</i> qPCR Assay | 96 Tests |
| QP10705 | HSV1/HSV2/TP <i>flex</i> qPCR Assay | 96 Tests |
| QP10706 | MG/TV <i>flex</i> qPCR Assay | 96 Tests |
| QP10707 | MH/UU <i>flex</i> qPCR Assay | 96 Tests |