

PRODUCT INFORMATION

EV4101, EV4102, EV4103, EV4104 and EV4105

Evidence and Evidence Investigator

The Kit Catalogue numbers above have undergone product improvement with the addition of the Benzodiazepines 3 assay to widen the detection of Benzodiazepine compounds, primarily Clonazepam and 7-Aminoclonazepam. In addition to this, these kits will be compatible with the newly recommended 1ng/ml cut off for Fentanyl.

These product improvements will be introduced from the below batch numbers onwards. Please deplete all current stock before proceeding with the batch numbers below, as previous lots of calibrators and controls will be incompatible with future batches of kit.

Please contact technical.services@randox.com if you require further information.

Device Name	Catalogue Number	Batch / Lot number
DOA Ultra Urine Evidence Full Kit	EV4101	4399 onwards
DOA Ultra Urine Evidence Half Kit	EV4102	4397 onwards
DOA Ultra Urine Investigator Kit	EV4103	4396 onwards
DOA Ultra Urine Calibrator Kit	EV4104	4163 onwards
DOA Ultra Urine Control Kit	EV4105	4400 onwards

PCA7054